Best in show
**AMERICAN ACADEMY
OF OPHTHALMOLOGY**
The Eye M.D. Association

D-EYE™ SMARTPHONE-BASED RETINAL SCREENING SYSTEM

D-EYE

INTRODUCING THE ONLY PORTABLE OPHTHALMOSCOPE CAPABLE OF RECORDING HIGH-DEFINITION IMAGES AND VIDEO OF THE RETINA

The innovative D-EYE Portable Retinal Imaging System easily attaches to a smartphone (Apple or Samsung), creating a ophthalmic camera capable of recording high-definition video and still images of the interior of the eye for health screening and evaluation.


A VALUABLE, MOBILE ACQUISITION TOOL

D-EYE can be used to conduct routine eye examinations and retinal screenings anywhere in order to allow clinician to:

- See the optic nerve head without dilating drops
- Look for neurological disorders associated with the eye
- Record images of children and infants without dilation drops
- Use to examine bed-ridden patients
- Use in emergency room evaluations
- Review before and after treatment images with patients, and parents or guardians

USED TO RECORD IMAGES AND VIDEO OF


GLAUCOMA

MACULAR HOLE

ASTROCYTIC
HAMARTOMA

OPTIC GLIOMA

UNDILATED
CHILD

D-EYE FEATURES

- Field of view of up to 20 degrees
- Portable - requires no external power or lighting
- Prescreen anywhere and share with other specialists immediately via the web
- Non-invasive, convenient and ergonomic
- Eliminates corneal glare
- Provides high-definition images
- Record multiple images or videos
- Built-in patient information file
- Optional private and secure cloud-based storage system
- iOS and Android based applications

CERTIFICATION AND PRIVACY

The D-EYE organization fulfills the Quality Management Requirements of ISO 13485, FDA GMP and HIPAA Patient Privacy Standards.


Required product safety and certifications will be completed where required for global acceptance in the medical community.

Current registrations: U.S. Food and Drug Administration - PJZ Ophthalmic camera regulation number 886.1120 class 2 510(K) exempt; CE - Medical Device Directive 93/42/CEE class I
D-EYE may be classified as an ISO 15004-2, group 2, device with some smartphones.

APP APPLE iOS AND ANDROID

D-EYE is easy to learn and simple to use. Download the smartphone app, affix the mounting bumper and fundus lens, set up patient/exam file, focus the phone's camera and begin recording. It's that easy, that simple.

Download the D-EYE app:

D•EYE

D-EYE App Guide

App Version	
iPhone	1.1.12
Android	1.1.0

Description

The portable D-EYE retinal and eye examiner attaches to a smartphone to record high-resolution images and video of the interior of the eye for ophthalmic and clinical assessment. The innovative fundus lens system allows regular screenings of the eye providing specific information about noticeable eye diseases, and it also can be used as a pre-screening recording device for further evaluation of specific medical conditions.

The low-cost D-EYE system will enable sharing the images and related data via the HIPAA-compliant D-EYE ImageVault™ cloud service (available late 2015) for diagnosis by various medical specialists, facilitating wide-spread screening anywhere, even in developing countries, bringing critical eye screening services to people in the most remote locations.

WARNINGS

- For safety reasons the D-EYE App does not work on jailbreak Apple smartphone.
- Install 1.1.9 and later versions for Apple devices with iOS 9.
- The supported smartphones are Samsung S4 and S5 and iPhone 5/5S and 6.

D-EYE Srl

Via Niccolò Tommaseo 77
35131 Padova, (PD), Italy
Tel. +39 049 781 1860
www.d-eyecare.com
info@d-eyecare.com

The final design, including color, may differ from printed images.

© Copyright 2015. D-EYE Srl. All Rights Reserved.

D-EYE, D-EYE Vet and the D-EYE Portable Retinal Imaging System are trademarks of D-EYE Srl.

SIGN IN AND LOGIN




D-EYE is easy to learn and simple to use. Download the smartphone app, affix the mounting bumper and fundus lens, set up patient/exam file, focus the phone's camera and begin recording. It's that easy, that simple. You can download the D-EYE app for your smartphone from either the Apple Store or Google Play Store. It can be found by searching for D-EYE. Once downloaded, you can launch the D-EYE app and register.

- 1 READ THE DISCLAIMER CAREFULLY. You must accept and confirm the rules and the restrictions of the disclaimer to access the D-EYE login page.
- 2 Log into the D-EYE app. If you have not registered, click on Signup in the upper left.
- 3 Provide a user name and password for your account and click Signup to get the D-EYE app. The password can be reset, but previously acquired data will be lost.

GLOSSARY


For health professionals who may not be familiar with all diseases and conditions of the eye, the app contains a pathology glossary, with images, to assist in evaluating the photos and video recorded by the D-EYE.

- 4 A pathology glossary is available in the app.
- 5 Select a pathology for more information.
- 6 Swipe left and right to scroll through pages.

RETINA SCREENING


Verify that the mounting bumper and D-EYE lens are firmly attached to the phone and start the examination. To begin, click on “Retinal Screening” and set up the parameters for the examination.

- 7 Select Retinal Screening.
- 8 Using the touch-screen keypad, enter the patient's name or an identification code.
- 9 Check the settings by clicking on “Settings” in the upper right corner of the screen, and change as needed.

Retina Screening - SETTINGS


The D-EYE smartphone app allows you to manage various settings for the acquisition of images and video.

Once you have adjusted the settings for the examination, click “Settings” in the upper left corner to return to the “Retinal Screening” page to begin the exam.


FOCUS PROCESS

This application requires some information of the patient in order to configure the focus optimally.

HOW TO CONFIGURE FOCUS

You can set the camera focus step-by-step before starting the examination, by following the guide on the screen of the smartphone. The app asks some questions about dilation and refractive error, and then indicates on how to set the focus. If you are an advanced user, you can use the manual configuration of the focus.

10 Enter the patient's name or ID code, configure the examination for either left or right eye viewing, select multi-shot or video mode and press "Start Examination" to begin the setup of the focus.

A Without error selected:

1. Attach the D-EYE device
2. Focus on a target at an infinite distance (>3 meters) and lock the focus

B With hyperopic error selected (pupils less than 3mm):

1. Attach the D-EYE device
2. Focus on a target at a distance according to the hyperopic error selected (example 12 inches/30cm at 3 diopters)
3. Lock the focus.
4. Perform the examination.




C With myopic error selected (pupils less than 3mm):

1. Remove the D-EYE lens from the phone
2. Focus on a target at a distance according to the myopic error (example 7 inches/17cm at 4 diopters)
3. Lock focus.
4. Reconnect the D-EYE lens.
5. Perform the examination.

11 For dilated eye, start the examination. For undilated eye, select and configure the degree of the refractive error.

Retina Screening - MULTISHOT



The application allows you to record the examination either as a series of still images (Multishot mode) or a full motion video. In Multishot mode, the camera will take a sequence of still shot images. You can control the number of images taken and the time between shots in the "Settings" section of the app. For Android device, there are not the sliders for the light because there's no need to adjust them.

- 12** After entering the patient's name or ID code, select "Multishot" mode and the eye you are about to examine, right or left.
- 13** Before performing the examination configure the settings of LED clicking on the  icon.
- 14** Click on the slider and the button near the  icon to adjust the intensity of the LED or to turn it off. Click on the slider near the  icon to adjust the exposure for the contrast.

Click on the gray circle to start recording and wait for the acquisition of a sequence of images. Adjustments of light and exposure are stored for the next examination.


If you choose the manual configuration of the focus, the  icon appears that allows you to set the focus manually.

- 15** When the Multishot photo sequence is complete, you will see the "Result" page and the multiple images that were taken during the exam. To view an enlarged image of a photo, press on the image until the larger version appears on the screen.
- 16** Pressing on the thumbnail image brings up the larger image on the "Analyze" page, allowing you to evaluate the image and delete any that are not suitable for diagnosis.

Retina Screening - MULTISHOT

Next, you can choose what to do with the examination result. The D-EYE app allows you to share your examination with the patient, other doctors and within your healthcare or hospital system. Images and video can be shared without Personal Health Information (PHI) via email or the phone's camera roll. Coming later in 2015, you will be able to upload images and video (with PHI) to the D-EYE ImageVault™, a secure HIPAA-compliant cloud service, to store the examination into the Patient Health Record.


A Restart the examination.

17 On the "Result" page, choose what to do with the examination result:


- A)** Clicking on "discard", restarts the examination.
- B)** Clicking on "save", archives the examination.
- C)** Clicking on the share icon  , the pictures will be sent by email, or saved in the phone's camera roll photo gallery.
- D)** Clicking on the file icon in the lower right of the screen  allows you to associate a pathology with the exam.

On the top it's possible to add some notes about the examination.

B The exam was saved in the archive. From here you can access the archive with your password or you can perform a new exam.


C You can easily share your examination with the patient, other doctors and within your healthcare or hospital system. With the D-EYE app, images and video can be shared (without PHI) via email or the phone's camera roll.


D Tag a pathology to the image(s). Associating a pathology with the exam helps to locate the exam results in the archive, and to classify the exam for discussion with other physicians.

Retina Screening - VIDEO


In addition to Multishot mode, the D-EYE app allows you to record the retina examination in full motion video. The video mode is helpful for obtaining a complete panoramic view of the posterior of the eye, as you are able to move the phone up and down, and side-to-side to capture a full view. Through post-exam editing, you can select the very best frames for archiving (and reducing the video file size). For Android device, there are not the sliders for the light because there's no need to adjust them.

18 Enter the patient's name or ID code, configure the examination for either left or right eye viewing, and select "Video" mode.


19 Before performing the examination configure the settings of LED clicking on the  icon.


20 Click on the slider and the button near the  icon to adjust the intensity of the LED or to turn it off. Click on the slider near the  icon to adjust the exposure for the contrast.

Quick Review of the Image

21 Before starting with recording, tap the screen for a quick review of the image.

22 With the still image you can zoom in with a double tap or a pinch. With a single tap you can close the quick review and come back to the examination.


Click on the red circle to start recording. The red circle will turn to a red square when the recording begins. To stop the recording, just press the red square. Adjustments of light and exposure are stored for the next examination.

If you choose the manual configuration of the focus, the  icon appears that allows you to set the focus manually.

ARCHIVE

Access the archive at any time to view saved examinations. Saving exams is useful in comparing images of similar pathologies and tracking the progression of a patient's disease or condition over time.


23 Examination data are encrypted in the smartphone. Without a password, the user cannot access the data. To view past, saved exams, click on "Archive" on the home screen.

24 Enter your password to access the archive.

25 Encrypted data of individual exams are sorted by pathology tag, date, patient name and note.

NOTIFICATIONS

The app is always in touch with you.


26 You can access the notifications available for you.

27 Read the notifications.

SUPPORT


Learn about D-EYE products, view online manuals, get the latest downloads, and more. Get service, support, and professional advice from the D-EYE team.

- 28** Clicking Support will open an email to customerservice@d-eyecare.com in which you can report your feedback or ask for assistance.
- 29** By selecting Support you open an email to us at customerservice@d-eyecare.com to send us:
- Feedback on reactions to the product (device + app);
 - Support, Technical Issues and troubleshooting requests;
 - or download the User manual.

SUPPORT / Assistance Mode

At the first launch of the D-EYE App, the assistance mode is active. It's possible to re-activate the assistance mode from the support menu.

- 30** Select "Assistant mode" on the support screen.

CAUTION

Read *“User Guide Annex I”* (p.153)

Leggi *“User Guide Annex I”*

Lire *“User Guide Annex I”*

Leer *“User Guide Annex I”*

Lesen *“User Guide Annex I”*

Ler *“User Guide Annex I”*

D•EYE

SDE1312W10001

English

Italiano

Français

Español

Deutsch

Português

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Pilipino

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D-EYE User Guide

**Errata corrigé de "D-EYE User
Guide REF:SDE1312W10001"**

*U.S. Food and Drug Administration
- PJZ Product Code substitutes
HKI on "Disclaimer" pages number:
7, 31, 43, 55, 67, 79, 91, 103, 115,
127, 139*

The final design, including color, may differ from printed images.

Copyright © 2015 D-EYE Srl. All Rights Reserved.
Printed in Italy

ACTIVATION CODE

•English	06
•Italiano	30
•Français	42
•Español	54
•Deutsch	66
•Português	78
•Русский	90
•中国	102
•日本人	114
•Pilipino	126
•한국의	138

• Disclaimer	07
• Mechanical Installation	08
• Light conditions for the exam	18
• Tips	19
• App Installation	20
• Focus	23
• Safety Information	24
• Security, Confidentiality	26
• Warranty	27
• Symbols List EN980	29

WARNING! USE OF THIS APP WITHOUT THE D-EYE INSTALLED IS STRICTLY PROHIBITED! THIS APPLICATION IS NOT INTENDED TO BE USED AS A STAND ALONE SOFTWARE. SERIOUS INJURY MAY OCCUR ALSO IF YOU USE A SMARTPHONE WITHOUT THE D-EYE DEVICE TOGETHER WITH APPLICATION. The product is to be used by a licensed healthcare practitioner or a licensed trained operator. The acquisition of the images must be always and only with the device hardware installed as instructed. Incorrect use of the application (no device connected or connected incorrectly) can cause damage to the device and to the people involved. Required product safety and certifications will be completed for global acceptance in the medical community. Current registrations:

- U.S. Food and Drug Administration : HKI Ophthalmic camera class II 510(K) exempt;
- Medical Device Directive 93/42/CEE class I.

D-EYE with some smartphones could be a ISO 15004-2 group 2 device.

CAUTION - The light emitted from this instrument is potentially hazardous. The longer the duration of exposure, the greater the risk of ocular damage. Exposure of light from this instrument when operated at maximum intensity will exceed the safety guideline after 3 min.

The range of focus adjustment for compensation of patient's refractive error is: -10 D to +8 D.

The D-EYE organization fulfils the Quality Management Requirements of ISO 13485, FDA GMP and HIPAA Patient Privacy Standards.


Complaints and reportable events may be provide to D-EYE by app/website section or by email to support@d-eyecare.com, or by mail to *D-EYE Srl via Niccolò Tommaseo 77 35131 Padua Italy*. Please contact us should you like to receive more information about the product and our current research programs.

EN

Upper side


Bottom side


1

Remove the D-EYE device from the box.


2

Take your smartphone bumper and open the screws or the clip.

3


Open the bumper and insert your smartphone, taking care to match the openings on the bumper with the battery connection and buttons of the phone.


4

Press the bumper on the smartphone and close the screws or the clip.

EN


iPhone 5/5S


iPhone 6

5

Place the bottom side of the D-EYE device in contact with the smartphone. A magnetic system will help correctly position the device on the smartphone. Look at the possible attached sheet to see the mechanical installation for the other smartphones.


Samsung Galaxy S4


Samsung Galaxy S5

EN


iPhone 5/5S


iPhone 6

6


The transparent circle on the bottom side of the D-EYE device should be aligned with the camera of the smartphone. The center of the transparent window should be aligned with the LED of the smartphone.

Samsung Galaxy S4


Samsung Galaxy S5

7

Start the D-EYE App. Position the smartphone near the patient's face (1-3 cm) with the D-EYE device facing the patient and the smartphone display facing the screener.

8

Start the examination of the *fundus oculi*.

EN

Exam is best conducted in dim light conditions. Panoramic retinal examinations (till the equator) can be performed only with dilated pupil. The optic nerve head and the posterior pole can be assessed through undilated pupil, with diameter not smaller than ~3 mm.

To assure that the patient's nose is not in the way of the camera, for iPhone models we suggest keeping the smartphone vertical for the right eye and horizontal for the left eye.


Acquisition is conducted similarly to direct ophthalmoscopy by placing the D-EYE at a distance of 1-3 cm from the patient's eye: the shorter the distance from the eye, the wider the observed field of view.

A video showing the correct acquisition procedure is available at: <http://www.d-eyecare.com>

For most eyes, suggested light intensity (option available for iPhone users only) is between 2/4 and 4/4 according to the size of the pupil and retinal pigmentation.


Tips: To visualize the optic nerve head, ask the patient to look straight ahead and approach with the phone slightly tilted toward the patient's nose. To visualize the macula, approach the eye straight on, toward the posterior pole, or ask the patient to look directly at the light from the D-EYE. To visualize the peripheral retina, tilt the smartphone toward the retinal quadrant to be examined and/or ask the patient to look in that direction (i.e. to examine the temporal retinal quadrant, slightly tilt the light beam temporally and ask the patient to look temporally).

We suggest selecting “video” retinal screening instead of “multishots” during the learning curve of the procedure.

9

Open the App Store and search for the D-EYE app. Install the application on your iPhone. Download the App Guide on <http://www.d-eyecare.com/documentation>.

10

Open the Play Store and search for the D-EYE app. Install the application on your Android phone. Download the App Guide on <http://www.d-eyecare.com/documentation>.

11

READ THE DISCLAIMER CAREFULLY. You must accept and confirm the rules and the restrictions of the disclaimer to access the login or signup D-EYE page.

- **Dilated eye:** point the smartphone at the patient's face and move the D-EYE near the eye, aiming at the red eye reflection. The autofocus will automatically compensate for the refractive error of the eye. Once the smartphone correctly focuses the retina, tap the  icon to lock in the correct focus. Perform the retinal examination.
- **Undilated eye with no refractive error:** for emmetropic undilated eyes. Keep the autofocus active, focus on a target at an infinite distance (~3 meters), Lock the focus with a tap on the  icon, Perform the examination.

▪ **Undilated eye (pupil's diameter smaller than 3 mm) with hyperopic error:**

1. Keep autofocus active.
2. Focus on a target at a distance according to the following table:

Hyperopic Error	Target Distance
2 diopters	50 cm / 20 inches
3 diopters	30 cm / 12 inches
4 diopters	25 cm / 10 inches
5 diopters	20 cm / 8 inches

3. Lock focus with a click on .
4. Perform examination.

▪ **Undilated eye (pupil's diameter smaller than 3 mm) with myopic error:**

1. Remove the D-EYE device.
2. Keep autofocus active.
3. Focus on a target at a distance according to the following table:

Myopic Error	Target Distance
2 diopters	12 cm / 5 inches
4 diopters	17 cm / 7 inches
6 diopters	25 cm / 10 inches
8 diopters	50 cm / 20 inches
10 diopters	infinite (~3 meters)

4. Lock focus with a click on .
5. Reconnect the D-EYE device.
6. Perform examination.

- **General Warning:** Failure to follow these instructions could result in injuries or damage to D-EYE or to people involved in using the D-EYE.
- **Handling and Transportation:** D-EYE is made of fragile components. Do not drop, disassemble, burn, or paint. Do not use if damaged. The upper and bottom face are made of glass. These parts can be damaged by a drop or a strong shock. Do not use a damaged D-EYE and do not touch the glass faces. Put D-EYE in its packaging after using.
- **Smartphone Compatibility:** D-EYE is designed and tested with compliant smartphones. Using the D-EYE with non-compliant smartphones could damage the D-EYE product, the smartphone and the people who are being screened; moreover D-EYE analysis results may not be valid or reliable.
- **Cleaning:** Clean the D-EYE metal ring and bumper with a soft cloth soaked in alcohol or mild cleaner/disinfecting solution. Clean the lenses with a clean optical grade cloth. Do not use any solvent-based cleaners on the lenses. Do not immerse D-EYE in any solution.
- **Disinfecting:** Do not sterilize the D-EYE.

- **Water:** Do not expose the D-EYE to water or rain. Avoid using and storing in moist places. Do not spill liquid or food on the D-EYE. Allow the D-EYE to air-dry; do not use heat sources.
- **Choking Hazard:** Some D-EYE accessories may present a choking hazard to small children. Keep these accessories away from small children.
- **Temperature Range:** The D-EYE is designed to operate from -20°C to +45°C. Operating D-EYE outside of this temperature range could damage the device. Avoid using the D-EYE during rapidly changing temperature or humidity conditions to prevent condensation over or inside lenses.
- **Magnets:** D-EYE uses 4 magnets to attach the device to the smartphone. Do not put magnets close to magnet-field-sensitive items such as credit and bank cards, pacemakers, defibrillators, electric and electronic systems, display screens, notebook, mechanical watches, or speakers. The magnet has a protective covering. Dropping or striking on hard surfaces could irreversibly damage the magnet. The covering contains nickel; avoid prolonged contact with skin.

1. Integrity, and Availability. Any password used for the D-EYE application is for individual use only. Users are responsible for the security of the password.
2. The security features of the D-EYE application enable users to meet obligations to ensure the confidentiality, integrity, and availability of patients' health information. Protecting patient data is solely the responsibility of the user. Users must verify patient identity before storing or sending patient records. D-EYE also recommends enabling a passcode in iOS to further protect the iPhone.
3. Managing data backup is the user's responsibility. D-EYE is not responsible for data backup, restoration or archiving.
4. Examination data are encrypted in the smartphone. Without a password, the user cannot access data. The password can be reset, but previously acquired data will be lost.

1. D-EYE Srl warrants the D-EYE to be free of original defects in material and workmanship and to perform in accordance with the manufacturer's specifications for a period of one year from the date of purchase. THIS WARRANTY APPLIES ONLY TO THE ORIGINAL CONSUMER PURCHASER OF THIS PRODUCT.
2. If this instrument is found to be defective or at variance with the manufacturer's specifications during the warranty period, D-EYE will repair or replace the instrument or component with a product of like kind or quality, which may be new or reconditioned with no charge to the purchaser for parts or labor. The repaired or replaced product will be warranted for 90 days from the date of return shipment or for the balance of the original warranty, whichever is longer.
3. This warranty only applies to D-EYE systems purchased new from D-EYE or its authorized distributors or representatives.
4. The purchaser must return the D-EYE system directly to D-EYE or an authorized distributor or representative and bear the cost of shipping. For repair or replacement on a D-EYE product the purchaser must obtain a Return Material Authorization (RMA) number from D-EYE support by submitting a repair request via email writing to support@d-eyecare.com.
5. This warranty does not cover breakage or failure due to tampering, misuse, neglect, accidents, modification or shipping.

EN

6. D-EYE SHALL NOT BE LIABLE TO YOU OR OTHER PERSONS FOR CLAIMS, LOSSES, COSTS, EXPENSES, DAMAGES WHATSOEVER DERIVING FROM MISUSE OF THE PRODUCT. D-EYE WILL, IN NO WAY, BE HELD RESPONSIBLE FOR ANY DAMAGES RESULTING FROM AN INCORRECT USE OF THE DEVICE EVEN IF D-EYE HAS WARNED ABOUT SUCH POSSIBILITIES.

7. This warranty is void if the instrument is not used in accordance with the manufacturer's recommendations or if repaired or serviced by anyone other than D-EYE or a D-EYE authorized representative. No other express or implied warranty is given.

Symbols	Title	Description
	CE marking of Conformity	Indicates the CE marking of Conformity as defined in EU directive 93/42/EEC, Annex XII
	Manufacturer	Indicates the medical device Manufacturer as defined in EU directive 93/42/EEC.
	Serial Number	This symbol shall be accompanied by the manufacturer's serial number.
	Catalogue Number	This symbol shall be accompanied by the catalogue number.
	Caution	Consult accompanying documents.

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- U.S. Food and Drug Administration: HKI Ophthalmic camera class II 510(K) exempt;
- Medical Device Directive 93/42/CEE class I.

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Copyright © 2015 D-EYE Srl . All Rights Reserved.

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Copyright © 2016 D-EYE Srl. All Rights Reserved.
Printed in Italy

D-EYE INTENDED USE

FOR THE ACQUISITION OF IMAGES OF THE EYE AND RETINA

MODELS AND COMPATIBLE SMARTPHONES

DEI0100W1S1400	iPhone 5/5S/SE
DEI0200W1S1400	iPhone 6/6S
DEI0300W1S1400	iPhone 6+/6S+
DES0100W1S1400	Samsung S4
DES0200W1S1400	Samsung S5


U.S. Food and Drug Administration - PJZ Ophthalmic camera
regulation number 886.1120 class 2 510(K) exempt.

D-EYE Srl Headquarter
Via Niccolò Tommaseo 77
35131 Padova, (PD), Italy
Tel. +39 049 781 1860
Fax. +39 049 781 1890
www.d-eyecare.com

EN

1

Place the bottom side of the D-EYE device in contact with the smartphone.
A magnetic system will help correctly position the device on the smartphone.

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Copyright © 2016 D-EYE Srl . All Rights Reserved.